# Anticorps monoclonaux anti- ***Enterocytozoon bienewsi*** et ***Encephalitozoon intestinalis***.

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## Références :

Simple species diagnosis of human intestinal microsporidia by an immunofluorescence test using specific monoclonal antibodies : an evaluation study in two hospitals in France. M.Thellier, I. Accoceberry, I. Desportes, S. Biligui, E. Bart-Delabesse, C. Ripert, M. Danis, A. Datry. First United Workshop on Microsporidia from Invertebrate and Vertebrate Hosts (NATO). July 12-15,2004. České Budějovice.

Evaluation of an immunofluorescent-antibody test using monoclonal antibodies directed against *Enterocytozoon bienewsi* and *Encephalitozoon intestinalis* for diagnosis of intestinal microsporidiosis in Bamako (Mali). Cisse O.A., Ouattara A., Thellier M., Accoceberry I., Biligui S., Minta D., Doumbo O., Desportes-Livage I., Thera M.A., Danis M., Datry A. The journal of Clinical Microbiology, 2002, **40** : 1715 -1718.

**BORDIER AFFINITY PRODUCTS SA**

Chatanerie 2, CH-1023 Crissier, Switzerland. Phone: + 41 21 633 31 67, <http://www.bordier.ch/>

# Monoclonal antibodies anti- *Enterocytozoon bienusi* and *Encephalitozoon intestinalis*

Technical sheet and instructions for article N° 8100, N° CE : H-CH/CA01/IVD/12173

## Intended use :

Routine stool examination for species specific diagnosis of **intestinal Microsporidia** by indirect immunofluorescent antibody assay (IFA) .

## Presentation :

The box contains 5 tubes:

**2 x 0.5 ml** of monoclonal antibody anti- *Enterocytozoon bienusi* ready to use (red caps).

**2 x 0.5 ml** of monoclonal antibody anti- *Encephalitozoon intestinalis* ready to use (green caps).

**1 x 2 ml** of ready to use Alexa fluor 488 labeled anti-mouse IgG conjugate containing Evans blue.

## Methode :

-Dilute feces with PBS (Phosphate-buffered Saline) (1 volume feces + 2 volumes PBS). Filter through 50 µm (best) or 100 µm filter.

-Put 2 µl of the fecal sample suspension to be tested on 18-well slides and dry for one hour.

-Fix the slides with methanol, let dry and dip them subsequently for 10 min in acetone at -20 °C., let dry.

-Rehydrate with a drop of PBS. Leave the slides for 5 min at room temperature, then aspirate the PBS on the side of the wells.

-Add 20 µl of either of the monoclonal antibodies on the wells and incubate for 30 min at room temperature in a humid atmosphere.

-After the 30 min incubation with the monoclonal antibodies wash 3 times with a drop of PBS, aspirating the PBS on the side of the wells.

-Add 20 µl of ready to use conjugate to each well. Incubate for 30 min at room temperature in a humid atmosphere, in the dark.

-Aspirate the conjugate, wash with a drop of PBS as specified above. Dip the slides in 3 changes of PBS.

-Drain the buffer from the slides. Dry the slides underneath and carefully around the samples, without touching the antigen.

-Add 3 drops of anti-fading fluorescence mounting medium, place a coverslip (24 x 60 mm) on the wells, avoiding the formation of air bubbles.

-Observe with a fluorescence microscope, equipped with the appropriate fluoresceine filter and an immersion objective (x 1000).The slides are stable several days at +4 °C in the dark.

## Interpretation :

To assess the specificity of the labeling, a positive and a negative sample as well as a sample without primary antibody (with the conjugate alone) are processed in parallel with the stool samples.

The monoclonal antibodies react exclusively with the spore walls of Microsporidia. *E. bienusi* spores (1.3 x 0.7 µm) and *E. intestinalis* spores (1.7 x 1.0 - 1.1 µm) are surface labeled with a marked peripheral fluorescence.

A sensibility close to 100% and a specificity better than those of specific staining methods (Weber and Uvitex 2B) were observed using the monoclonal antibodies and IFA.

## References:

Simple species diagnosis of human intestinal microsporidia by an immunofluorescence test using specific monoclonal antibodies : an evaluation study in two hospitals in France. M. Thellier, I. Accoceberry , I. Desportes, S. Biligui, E. Bart-Delabesse, C. Ripert, M. Danis, A. Datry. . First United Workshop on Microsporidia from Invertebrate and Vertebrate Hosts (NATO). July 12-15,2004. Ceské Budějovice.

Evaluation of an immunofluorescent-antibody test using monoclonal antibodies directed against *Enterocytozoon bienusi* and *Encephalitozoon intestinalis* for diagnosis of intestinal microsporidiosis in Bamako (Mali). Cisse O.A., Ouattara A., Thellier M., Accoceberry I., Biligui S., Minta D., Doumbo O., Desportes-Livage I., Thera M.A., Danis M., Datry A. The journal of Clinical Microbiology, 2002, **40** : 1715 - 1718.

## BORDIER AFFINITY PRODUCTS SA

Chatanerie 2, CH-1023 Crissier, Switzerland. Phone: + 41 21 633 31 67, <http://www.bordier.ch/>