

# Chemistry Model-based Method for Estimating DNA Amounts from PCR Amplified Samples

Benedict G Archer, Jules Chen, Mark Dobbs & Maurice Exner

Focus Diagnostics, R&D Molecular Products

## Background & Objective

Methods in common use for estimating target NA amounts in PCR-processed samples ignore the chemistry occurring each cycle. Rather PCR data is treated as resulting from a growth process, relying either on assumptions related to replication efficiency in cycles up to the first exceeding a (usually low) threshold signal, or on estimated parameters of a fitted logistic function. Our objective: to model the chemical reactions of template replication through enough cycles *beyond the geometric phase* of the PCR process to enable accurately estimating target NA amounts without resorting to special phenomenological models. Our starting point is reference [1].

## Extension Reactions


primer  $p$  binding to single strand template  $s$



polymerase  $q$  binding to primer-template duplex  $h_0$



primer extension in replication complex  $r_0$



amplicon-primer competition for template



## Approximations & Model Variables

Idealizing approximations:

- all annealed primers extend completely  $K_3 \approx \infty$
- primer in large excess  $[p] \approx [p]_{total} \quad K_1[p]_{total} \gg 1$

Model variables:

- reduced template concentration  $S = K_2[s]_T$
- reduced polymerase concentration  $Q = K_2[q]_T$
- amplicon-primer competition  $D = K_4/[K_2(K_1[p]_T)^2]$

## Model Parameters and Dependencies

parameter	meaning	computation mode			dependence
		full calibration	system calibration	prediction	
$X_0$	sample target NA amount	fixed	fixed	estimated	sample
$\mathcal{K}$	normalized $K_2$	estimated	fixed	fixed	chemistry
$Q$	reduced polymerase conc.	estimated	estimated	estimated	chemistry & sample
$D$	amplicon-primer competition	estimated	fixed	fixed	chemistry
$G$	unit conversion factor	estimated	estimated	fixed	system

## Model Summary

The measured net profile signal  $Y_k$  at cycle  $k$  is represented by a general phenomenological model including initial target amount  $X_0$ , a system gain parameter  $G$  and a product involving cycle efficiencies  $E_i$ ,

$$Y_k = GX_0 \prod_{i=1}^k (1 + E_i)$$

Within each cycle the PCR chemistry model is embodied in a polynomial in replication efficiency  $E_k = [r]_k/[s]_{T,k}$ , with coefficients functions of model parameters,

$$S_k^2 E_k^3 - (S_k^2 + (1 - D + 2Q)S_k)E_k^2 + Q(2S_k + Q + 1)E_k - Q^2 = 0.$$

A recursion links the reduced amplicon concentrations in the cycle sequence,

$$S_{k+1} = S_k(1 + E_k).$$

$S_0$  includes sample target amount  $X_0$  and  $\mathcal{K}$ ;  $\mathcal{N}$  and  $V$  denote Avogadro's number and reaction volume,

$$S_0 = X_0 \mathcal{K}. \quad \mathcal{K} = K_2/(\mathcal{N}V)$$

## Implementation

- amplification detection and NLLS background estimation
- select cycles from profiles having efficiency greater than defined fraction of maximum efficiency (Figure 1)
- assemble data for estimation of parameters common to all profiles by NLLS (Figure 2)




Figure 1. Selecting 10 cycles preceding first with below minimum efficiency




Figure 2. Profile segments from 3 replicates of 5 calibrators

## Conclusions

- A real model-based analysis of PCR data offers:
  - accurate, assumption-free, objective estimates of target NA amount
  - estimates of target NA amount uncertainty
  - precision as good or better than other methods
  - parameters having real, physical significance
  - reliable estimate of first cycle efficiency
  - useful validity checks
  - easy calibration transfer between systems

## References

- [1] D. Stolovitzky & G. Cecchi, PNAS **93** (1996), 12947-12952.

## Acknowledgements

bga thanks Ming-Chou Lee for encouraging this work and Jeffrey Adachi for insights and implementation advice and guidance.



## Comparison of chemistry model method with threshold method \*




Figure 3. Chemistry model analysis of CMV validation plot




Figure 4. Threshold method analysis of CMV calibration plot

CMV true $X_0$	chemistry model				threshold method			
	est $X_0$	%error	%CV	$E_1$	est $X_0$	%error	%CV	efficiency
200000000	199841706	-0.1	1.2	0.926	278136085	39.1	3.5	0.961
200000000	19331065	-3.3	1.0	0.925	24441198	22.2	4.6	0.961
2000000	1503930	-24.8	46.2	0.905	1598590	-20.1	59.2	0.961
200000	186896	-6.6	6.0	0.921	198072	-1.0	8.3	0.961
20000	20121	0.6	2.8	0.923	20541	2.7	3.8	0.961
2000	2053	2.7	6.6	0.923	1922	-3.9	5.3	0.961
200	220	10.1	24.2	0.921	186	-7.2	24.0	0.961
20	21	5.9	27.2	0.912	14	-29.2	36.7	0.961
2	14	587.5	8.1	0.896	6	210.5	2.7	0.961




Figure 5. Chemistry model analysis of EBV validation plot




Figure 6. Threshold method analysis of EBV calibration plot

EBV true $X_0$	chemistry model				threshold method			
	est $X_0$	%error	%CV	$E_1$	est $X_0$	%error	%CV	efficiency
200000000	208863846	4.4	23.5	0.916	240491477	20.2	16.1	0.973
200000000	18837207	-5.8	19.1	0.925	23109357	15.5	16.3	0.973
2000000	1836084	-8.2	22.5	0.926	1915625	-4.2	13.0	0.973
200000	184923	-7.5	9.5	0.925	187195	-6.4	8.2	0.973
20000	19295	-3.5	15.0	0.925	18295	-8.5	9.3	0.973
2000	2014	0.7	9.4	0.923	1678	-16.1	3.0	0.973
200	204	2.0	8.1	0.927	173	-13.6	6.0	0.973
20	21	4.9	8.1	0.922	15	-27.2	11.6	0.973
2	5	162.4	5.0	0.926	4	93.7	10.7	0.973




Figure 7. Chemistry model analysis of FLU validation plot




Figure 8. Threshold method analysis of FLU calibration plot

FLU true $X_0$	chemistry model				threshold method			
	est $X_0$	%error	%CV	$E_1$	est $X_0$	%error	%CV	efficiency
200000000	219835414	9.9	14.1	0.926	273152958	36.6	15.0	0.984
200000000	18158971	-9.2	14.5	0.922	18816832	-5.9	19.2	0.984
2000000	1638426	-18.1	9.3	0.930	1789355	-10.5	12.4	0.984
200000	157962	-21.0	6.5	0.931	158050	-21.0	9.5	0.984
20000	22789	14.0	20.1	0.931	20050	0.2	11.9	0.984
2000	2347	17.4	10.4	0.931	1967	-1.6	10.8	0.984
200	211	5.7	18.8	0.936	182	-9.2	17.3	0.984
20	27	33.2	41.0	0.945	27	33.2	30.8	0.984
2	2.4	20.0	—	0.928	n.d.	—	—	—
0.2	0.9	361.0	111.0	0.883	n.d.	—	—	—

\* samples processed on the 3M Integrated Cycler®